

Технологии медицинских имплантов

Фирма	Система имплантов	Платформа	Усилие	Тип посадки	Pre-mill	Scan	Ti-База	Аналог		
Straumann (ITI)	Bone Level	3.5 (NC)	35 Ncm	Engaging	4035-90B	4035-80BN	4035-70B	4035-40B		
				Non-engaging	4035-90BN		4035-70BN			
		4.1 (RC)		Engaging	4043-90B	4043-80B	4043-70B	4043-40B		
				Non-engaging	4043-90BN		4043-70BN			
	Tissue Level	4.8 (RN)		Engaging	4048-90	4048-80	4048-70	4048-40		
				Non-engaging			4048-70N			
		6.5 (WN)		Engaging	4065-90	4065-80	4065-70	4065-40		
				Non-engaging			4065-70N			
Nobel Biocare	Active	3.5 (NP)	35 Ncm	Engaging	5035-90A	5035-80A	5035-70A	5035-40A		
		4.3 (RP)		Engaging	5043-90A	5043-80A	5043-70A	5043-40A		
		5.0 (WP)		Engaging		5050-80A	5050-70A	5050-40A		
	Branemark	3.3 (NP)		Engaging	3035-90B	3035-80B				
		3.75/4.0 (RP)		Engaging	3041-90B	3041-80B	3041-70B			
		5.0 (WP)		Engaging	3050-90B	3050-80B				
	Replace	3.5 (NP)		Engaging	5035-90R	5035-80R	5035-70R	5035-40R		
		3.5 (NP)		Non-engaging			5035-70RN			
		4.3 (RP)		Engaging	5043-90R	5043-80R	5043-70R	5043-40R		
		4.3 (RP)		Non-engaging			5043-70RN			
		5.0 (WP)		Engaging	5050-90R	5050-80R	5050-70R	5050-40R		
		5.0 (WP)		Non-engaging			5050-70RN			
		6.0		Engaging	5060-90R	5060-80R	5060-70R	5060-40R		
	6.0	Non-engaging								
	Dentsply	Ankylos		Balans C/X	15 Ncm	Non-engaging	1100-90AN	1100-80AN	1100-60AN	
		Xive		3.0	24 Ncm	Engaging	7030-90	7030-80	7030-70	
3.4			Engaging	7034-90		7034-80	7034-70			
3.8			Engaging	7038-90		7038-80	7038-70			
4.5			Engaging	7045-90		7045-80	7045-70			
5.4			Engaging	7054-90		7054-80	7054-70			
Biomet 3i	Certain	3.25 (3.4)	20 Ncm	Engaging	6034-90	6034-80	6034-70	6034-40		
		4.0		Engaging	6041-90	6041-80	6041-70	6041-40		
		5.0		Engaging	6050-90	6050-80	6050-70	6050-40		
		6.0		Engaging	6060-90	6060-80				
	Osseo Tite	3.5		Engaging		3035-80T	3035-70T			
		4.1		Engaging	3041-90B	3041-80B	3041-70B			
		5.0		Engaging	3050-90T	3050-80T				
				Non-engaging						
Osstem	GS/TS	M	30 Ncm	Engaging	1100-90OM	1100-80OM	1100-70OM	1100-40OM		
				Non-engaging	1100-90OMN		1100-70OMN			
		R		Engaging	1100-90OR	1100-80OR	1100-70OR	1100-40OR		
				Non-engaging	1100-90ORN		1100-70ORN			
	SS	4.8		Engaging	4048-90OR	4048-80OR	4048-70OR			
		6.0		Engaging	4060-90OW	4060-80OW	4060-70OW			
Dentium	Superline Implantium	R	20-30 Ncm	Engaging	1100-90DR	1100-80DR	1100-70DR	1100-40DR		
				Non-engaging	1100-90DRN		1100-70DRN			
Dio	SM	M	30 Ncm	Engaging	1100-90TM	1100-80TM	1100-70TM			
				Non-engaging						
		R		Engaging	1100-90TR	1100-80TR	1100-70TR			
				Non-engaging						
Astra Tech	Osseo Speed	Yellow 3.0	25 Ncm	Engaging	1030-90	1030-80	1030-70			
		Aqua 3.5		Engaging	1035-90	1035-80	1035-70	1035-40		
		Lilac 4.5		Engaging	1045-90	1045-80	1045-70	1045-40		
Zimmer	Ad-Vent			Engaging	2045-90A	2045-80A				
	TSV	3.5	20-30 Ncm	Engaging	2035-90	2035-80	2035-70	2035-40		
		3.5		Non-engaging			2035-70N			
		4.5	Engaging	2045-90	2045-80	2045-70	2045-40			
			Non-engaging			2045-70N				
		5.7	Engaging	2057-90	2057-80	2057-70	2057-40			
			Non-engaging							
	Spline	3.75/4.0			Engaging					
5.0				Non-engaging						

Фирма	Система имплантов	Платформа	Усилие	Тип посадки	Pre-mill	Scan	Ti-База	ССМ - литье
Camlog	Screw Line	3.8	20 Ncm	Engaging	5038-90C	5038-80C		
		4.3		Engaging	5043-90C	5043-80C		
		5.0		Engaging	5050-90C	5050-80C		
Thommen Medical	SPI	3.5	25 Ncm	Engaging	6035-90S	6035-80S	6035-70S	
		4.0		Engaging	6040-90S	6040-80S	6040-70S	
		4.5		Engaging	6045-90S	6045-80S	6045-70S	
		5.0		Engaging	6050-90S	6050-80S	6050-70S	
		6.0		Engaging	6060-90S	6060-80S	6060-70S	
Bicon	Narrow	2.0 (M)		Engaging	1000-90M	1000-80M		
	Max	2.5 (R)		Engaging	1000-90R	1000-80R		
	Short	3.0 (W)		Engaging	1000-90W	1000-80W		
Kyocera	POI	3.7	20 Ncm	Engaging	7037-90J	7037-80J		
		4.2		Engaging	7042-90J	7042-80J		
		4.7		Engaging	7047-90J	7047-80J		
		5.2		Engaging	7052-90J	7052-80J		
BBC	Mytis Arrow	B Type	30 Ncm	Engaging	4048-90MYT	4048-80MYT		
		C Type		Engaging	1100-90MYT	1100-80MYT		
		E Type		Engaging	4043-90MYT	4043-80MYT		
Megagen	EZ Plus	M	35 Ncm	Engaging	1100-90MM	1100-80MM	1100-70MM	
		R		Engaging	1100-90MR	1100-80MR	1100-70MR	
		W		Engaging	1100-90MW	1100-80MW	1100-70MW	
	Anyridge	R		Engaging	1100-90MEA	1100-80MEA	1100-70MEA	1100-40MEA
		Non-engaging			1100-70MEAN			
	Anyone	R		Engaging	1100-90MAO	1100-80MAO	1100-70MAO	
Rescure	R	45 Ncm	Engaging	3050-90MRE	3050-80MRE			
Shinhung	Luna	R	30 Ncm	Engaging	1100-90SMR	1100-80SMR		
NeoBiotech	CMI	R	30 Ncm	Engaging	1100-90NER	1100-80NER	1100-70NER	
Dentis	One Q	R	30 Ncm	Engaging	1100-90DSR	1100-80DSR	1100-70DSR	
Warantec	Oneplant	R	30 Ncm	Engaging	4030-90WR	4030-80WR		
MIS	SEVEN	Narrow	35 Ncm	Engaging	E2032-90M	E2032-80M	E2032-70M	E2032-40M
		Standart		Engaging	E2037-90M	E2037-80M	E2037-70M	
		Wide		Engaging	E2044-90M	E2044-80M	E2044-70M	
	C1	Narrow		Engaging	E8027-90M	E8027-80M	E8027-70M	
		Standart		Engaging	E8032-90M	E8032-80M	E8032-70M	E8032-40M
		Wide		Engaging	E8040-90M	E8040-80M	E8040-70M	E8040-40M
SIC	ace 3.4 / 4.0 & max 3.7 / 4.2		30 Ncm	Engaging	6033-90SI	6033-80SI	6033-70SI	
	ace 4.5 / 4.0 & max 4.7 / 5.2			Engaging	6042-90SI	6042-80SI	6042-70SI	
BEGO	Semados S/ RI/ Mini Line	3.74	30 Ncm	Engaging	2037-90B	2037-80B	2037-70B	
		4.1		Engaging	2041-90B	2041-80B	2041-70B	
		4.5		Engaging	2045-90B	2045-80B	2045-70B	
		5.5		Engaging	2055-90B	2055-80B	2055-70B	
Anthogyr	Axiom		25 Ncm	Engaging	1100-90AT	1100-80AT	1100-70AT	
ICX	Templant			Engaging	1100-90ICX	1100-80ICX	1100-70ICX	

Тел. 8-499-714-70-30

Тел. 8-903-225-49-70