













HATHOfex S plus P HATHOfex S plus P



Schleif- und Polierwerkzeuge aus elastisch gebundenem Korund
Grinding and Polishing Instruments with elastic bonded corundum
Meule à affiler et à polir avec corindon, lien élastique
Herramientas de pulido de corindón elástico prensado
Abrasivi e strumenti per lucidare in lega di corindone

1:1

Form shape	Körnung grit	80	180	240	400	Ø mm	H mm	
		200 6 HP	210 6 HP	220 6 HP	230 6 HP	20	3	
∪ 20.000								
		205 17 HP	215 17 HP	225 17 HP	235 17 HP	20	3	
∪ 5.000								
		204 20 HP	214 20 HP	224 20 HP	234 20 HP	20	3	
∪ 15.000								
		202 6 HP	212 6 HP	222 6 HP	232 6 HP	6	16	
∪ 20.000								
		202 8 HP	212 8 HP	222 8 HP	232 8 HP	8	16	
∪ 15.000								