



T&F Hybrid Points

Ultrafine diamonds for finishing of tooth enamel, porcelain and composite

T&F Hybrid Points with their ultrafine diamond particles fulfill the requirements of modern, aesthetic dentistry regarding minimal corrections and best preparation for polishing.

The instruments are available in ten different sizes and shapes for the different preparation steps and are individually packed in plastic containers.

PN	0925-1	0927-1	0929-1	0931-1	0933-1	0935-1	0937-1	0939-1	0941-1	0943-1
ISO No.	018	014	018	014	018	012	009	012	013	011
Head length (mm)	1.6	3.0	3.5	3.0	3.5	4.0	3.7	4.2	6.0	8.0
Bur No.	7006	7104	7106	7404	7406	7803	7901	7903	LT1	LT2


Indications

CHAIRSIDE

Finishing of composites; very fine adjustments to tooth enamel and porcelain

Range of Application

Suitable for the following materials:
Tooth enamel, porcelain, composites



The T&F Hybrid Points Kit contains all 10 shapes.



Min⁻¹ Recom. 300,000 - 450,000 / CP 0 - 0.8 N

PN	Quantity	Product
0925-1 - 0943-1	6 pcs.	T&F Hybrid Points, FG shank
0950	1 pc.	T&F Hybrid Points Kit

BENEFITS

Exact, even distribution of very fine diamond particles

Computerised manufacturing technique

Consistent performance

Ideal surface finish before final polishing

Precise concentricity